

The 1st Chair Of Dentist's Clinic

Hot Sale Model:A001



◆ The Appropriate Functional Design



Aluminum Strong and weak suction

Illuminance degree 8000– 20000LUX
Color temperature 4500 –5000K – Yellow :3000K
Light spot size 80X160mm (distance:700mm)
Control Method: Induction Switch, Manual Control



Comfortable enlarged dental chair base and backrest ,bring safety and comfort to patients while lying on the dental chair



Three fold assistant tray , save the space of clinic room



Highly instrument tray with 6 holes

Assistant stool with metal base

The multifunctional foot pedal can control separately the water and air

Can control the dental chair and the backrest up and down



Adjustable armrests to meet operational needs



Big tray accommodates as much as surgical tools.

◆ Standard Configuration

Flash and adjustable lamp arm

6 holes LED sensor lamp

Microfiber Leather

Big 6 holes down-tray

Durable water and air tubing

Removable ceramic spittoon

Aluminum strong and weak suction head

Ceramic Removable Spittoon

Double armrest

Assistant stool with metal legs

Folding assistant frame

Four-function pedal

Weight : 200 KGS

Size : 1.00*1.30*0.85m /1.00*1.30*0.75m (without wooden legs)

Container max quantity :

20GP=20PCS / 20GP= 27PCS (without wooden legs)

40HQ=63PCS

◆ Details in the details



User-centric Design

Big tray with 6 holes, accommodates as much as surgical tools



Smart Dental Light

6 Holes , Induction Switch, Manual Control



Durable assistant stool

Assistant stool with metal legs , more valuable and stable



3 fold assistant tray Design

save the space of clinic room

◆ Standard configuration

PU Color for reference


ASK001
(A001-R)

ASK002
(A001-B)

ASK003
(A001-G)

ASK004
(A001-D)

◆ Installation Drawing



EC REP

Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175 Fax: +49-40-255726

Manufacturer: Guangdong JINME Medical Technology Co., Ltd.
Tel: +86-757-88777991
www.appldentalmed.com

2023/08 Version



◆ Additional Configuration



◆ One-Station Purchasing Options (Agents Fulfill the Additional Set Retail Price)

L5 +



L4 +



L3 +



L2 +



L1 +



◆ Standard Configuration and Optional Equipment



Appledental curing light

1. Dimensions: 26mm×29mm×240mm
2. Net weight: 120g
3. Power supply: 24V-50Hz/60Hz
4. Applied part: optical fiber
5. Light source: Blue light, Wave length: 420nm-480nm
6. Working condition: Environment temperature: 5°C-40°C Relative humidity: <80% Atmosphere pressure: 70kPa-105kPa
7. Consumption power: <8W
8. Protection type against electrical shock: class II
9. Protection against electrical shock: type B
10. Protection against harmful ingress of water or particular matter: ordinary equipment (IPX0)
11. Safety in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide: not suitable under this condition.



LED scaler

Compatible Brand: EMS
 Handpiece: HP-5L(with LED lamp, Detachable)
 Function: G, P
 Potentiometer: Carbon-film potentiometer (circular)
 Size (mm): 75*56*34mm
 Weight of main unit: 0.125kg
 Power input: 24VAC 50Hz/60Hz1.3A
 Fuse of main unit: T1.6AL 250V
 Primary tip vibration excursion: <200µm
 Tip vibration frequency: 28kHz±3kHz
 Output power of tip: 3W-20W
 Half-excursion force: 0.5N-5N
 Water entrance pressure: 0.01-0.5MPa



1 drive 1 compressor

MODEL: AP800
 Machine dimension: 400*400*630mm
 Rated Power: 560W
 Nominal: 70L/min
 Totation: 1380r/min
 Capacity: 32L
 Noise: 65db
 Weight: 30.5kg
 Carton Size: 465mm*465mm*695mm



WIFI Intra oral camera

Gross Weight: 6.0kg
 Monitor: 38*37*2cm
 LCD Clamp: 34*27*2cm
 Unit Package Size: 571cm*41.9cm*11.8cm
 Focal length: 3mm-50mm, Auto focus
Features:
 1. Multifunctional-Can take 1&4 pictures
 2. Clearer-High resolution SONY CCD 12 mega pixel
 3. Durable quality- HUA WEI Hisilicon chips
 4. Easy connections-Fast-plug connectors
 5. 8G storage-In USB disk
 6. HDMI+USB interface;
 7. 6-LEDS(with glass cover)
 8. Wifi connection



2pcs A1+1set low speed (CA+SH+AM)

Code: GREEN-L/BLUE-L
 Air Pressure: 0.25(Mpa)-0.3(Mpa)
 Water pressure: 198kpa(2kg)
 Forward rotation speed: 18 (krpm)
 Reversal rotation speed: 16(krpm)
 Torque: 0.4(N) -0.45 (N)
 Bursize: Ø2.334-Ø2.355(mm)
 Sterilization temperature: 135°C
 Waterflow: 90(min/ml)-110(min/ml)
 Water pressure: 0.2(Mpa)-0.23(Mpa)
 Noise: <65dB
 Chuck type: latch type
 Description: external water course, latch type, can use polisher, burs, dentallab disc, reamer, self lubricated
 function of air motor: 4 holes and 2 holes optional



Whitening machine

- (1). Power supply: AC100-240V/50/60HZ
- (2). Power: 300W
- (3). Clamping size: 45-50cm.
- (4). Light size: 18x74mm2
- (5). Light intensity: >6000mw/cm
- (6). Light timing: 0-20mins
- (7). Power instability: <1%
- (8). Working method: time output & continuous working
- (9). Coupling: pd yester monowave
- (10). Gross weight: 1.4kg/pcs